



**HA TECH**  
BioPharma & Diagnostics

# RESPIRATORY COMBO PANEL SARS-CoV-2

- COVID & Influenza A/B Rapid Antigen Test Kit
- Respiratory Syncytial Virus RSV (Self-Test)

**HA Tech introduces a respiratory panel which aids the diagnosis of 4 diseases at a time.**

## PRODUCT PURPOSE

SARS-CoV-2/ Inf A/ Inf B/ RSV Respiratory Self-Test Panel (RAT-19) is intended to aid the diagnosis of COVID-19, Influenza A, Influenza B or Respiratory Syncytial Virus (RSV) in symptomatic patients. This panel is designed to be tested by using just one nasal sample.

## FEATURES



High Sensitivity and Specificity  
•96% and 100% respectively  
Detect all variants



LoD of 40 TCID<sub>50</sub>/mL for SARS-CoV-2, 1.5 X 10<sup>3</sup> TCID<sub>50</sub>/mL for Influenza A, 5.0 X 10<sup>3</sup> Influenza B. 159 TCID<sub>50</sub>/mL for RSV



Simple Self-Testing (Nasal swab)




Fast Results & Easy Interpretation



Validated for extreme local weather conditions



Made in Australia



ORDER CODE

DESCRIPTION

QTY/Box

HA-ICR19P5

SARS-CoV-2 / Influenza A / Influenza B / RSV  
Respiratory Self-Test Panel

5 Combo  
Tests



+61 02 867 80993  
sales@wishmed.com.au  
www.wishmed.com.au