

Disposable or Re-usable INOCULATION LOOPS

RIGID or FLEXIBLE LOOPS

CALIBRATED PLASTIC LOOPS

**NICHROME NON-CALIBRATED
WIRE METAL LOOPS**



**For general
microbiology
and laboratories**



DISPOSABLE PLASTIC LOOPS HARD or SOFT/FLEXIBLE

- Diaguru Plastic inoculating loops provide regular wetting and full liquid transfer.
- Available with rounded tips.
- Each inoculating loop has smooth edges.
- Colour coded.

Features

Smooth and flexible for uniform, smooth streaking without damaging the gel surface.

Rounded tip.

Sterilized by Ethylene Oxide.

Sterile & disposable to eliminate the risk of cross-contamination & free from lubricants - does not contaminate any media.

HARD



1µl
D-INL1HBX

SOFT/FLEXIBLE




1µl
D-INL1SBX



10µl
D-INL10HBX



10µl
D-INL10SBX



METAL INOCULATION LOOP FOR NEEDLE HOLDER

MF6624000-EE

- Made of stainless steel wire.
- Length: approx. 35mm.
- Diameter of loop: approx. 5mm.
- Diameter of wire: approx. 0.5mm.
- 100pcs/Poly Bag.
- Capacity of loop approx. 7ul.

RE USABLE WIRE METAL LOOPS

- Diaguru Platinum inoculating loops are medical devices used for bacteriological and other microbiological applications.
- They will better resist constant exposure to acid and flame than the nichrome wire inoculating loop.
- Platinum Diaguru, wire inoculating loops, are faster cooling, very little waiting time.
- Resistance to oxidation and degradation from chemicals.
- Durable and long lasting.
- Calibrated and come with a calibration certificate.



2.5µl
D-INLM25EE



4.0µl
D-INLM40EE



Features

Designed and moulded with precision for accurate volumetric technique every time.

Validated to ISO standards.

Ensure better wetting with precise moulded smooth edges.

RE USABLE HANDLE for METAL LOOPS

Handle Features

The inoculating loop or needles are handled by holding the handle part of the loop or needle and facilitating fatigue-free application.

The inoculation holder is made of bakelite aluminum handle and brass rod in either an insulated or non-insulated format.

The handle is 8 inches long.

Shaft: Nickel-plated brass is used to make the loop shaft.

D-INLHANE



CODE	PRODUCT DESCRIPTION	PACK SIZE
Inoculation Loop + Incorporated Needle - Hard		
D-INL1HBX	Inoculation Loops, 1 μ L + Incorporated Needle, AS Material, Neutral Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile.	1000 pcs/box
D-INL1HCTN	Inoculation Loops, 1 μ L + Incorporated Needle, AS Material, Neutral Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile.	10 boxes/Ctn
Inoculation Loop + Incorporated Needle - Flexible		
D-INL1SBX	Inoculation Loops, Flexible 1 μ L + Incorporated Needle, ABS Material, Neutral Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile.	1000 pcs/box
D-INL1SCTN	Inoculation Loops, Flexible 1 μ L + Incorporated Needle, ABS Material, Neutral Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile	10 boxes/Ctn
Inoculation Loop + Incorporated Needle - Hard		
D-INL10HBX	Inoculation Loops, Hard 1 μ L + Incorporated Needle, AS Material, Blue Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile.	1000 pcs/box
D-INL10HCTN	Inoculation Loops, Hard 1 μ L + Incorporated Needle, AS Material, Blue Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile	10 boxes/Ctn
Inoculation Loop + Incorporated Needle - Flexible		
D-INL10SBX	Inoculation Loops, Flexible 10 μ L + Incorporated Needle, ABS Material, Blue Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile.	1000 pcs/box
D-INL10SCTN	Inoculation Loops, Flexible 10 μ L + Incorporated Needle, ABS Material, Blue Colour 20pcs/Zip-Lock Bag, 50pks/box, Sterile	10 boxes/Ctn
Inoculation Loop + Incorporated Needle - Metal		
D-INLM25EE	Inoculation Loops, Wire Loop End, Nickel Alloy with Chromium Material, Dia 2.5mm, Overall Length 50mm each.	1000 pcs/box
D-INLM40EE	Inoculation Loops, Wire Loop End, Nickel Alloy with Chromium Material, Dia 4.0mm, Overall Length 50mm each.	10 boxes/Ctn
MF6624000-EE	Marienfeld Inoculation Loops Metal, For Needle Holder, Dia 5.0mm, 7ul, 35mm length thick.	100 pcs/bag
Plastic Handle + 2 Inoculation Loops		
D-INLHANE	Each Plastic Handle, Overall Length 200mm, Supplied with Two Circle Ends Loops made of Nickel Alloy Material.	Combo/bag

PLEASE CONTACT US TO PLACE YOUR ORDER:

Phone +61 2 86786499

sales@diaguru.com.au

www.wishmed.com.au

