

Diaguru Orbital Shaker (D-SHAKE)

Code: DOS20P



- Platform shaker with true Orbital and Linear reciprocating motion for shaking mass upto 7.5 kg
- Powerful BLDC Motor and Microprocessor Controlled for Smooth operations and constant speed.
- Timer ranging from 1 min to 99 hours and 59 minutes
- User Friendly and Programmable: User can save upto 6 programs with each program having 9 sub programs, Set motion direction (Clockwise and Counter clockwise), speed, time dwell time, sequence and loop operation
- Different platform attachments makes it possible to use almost all shapes and sizes of vessels
- Option for auto power On (in case of power failure) – Resumes operation for the remaining time.

Optional Attachment



Universal




Funnel



Plastic Holder



SS Holder



Plate

ORDERING INFORMATION

D-Shake DOS10P / 20

Product Name	Cat #	Details
DOS10P	DOSP-1000	Platform Shaker with 10 mm orbital motion
DOS10P	DOSP-1001	Platform Shaker with 10 mm orbital motion along with RS 232
DOS10P	DOSP-1010	Platform Shaker with 10 mm linear motion
DOS10P	DOSP-1011	Platform Shaker with 10 mm linear motion along with RS 232
DOS20P	DOSP -2000	Platform Shaker with 20 mm orbital motion
DOS20P	DOSP -2001	Platform Shaker with 20 mm orbital motion along with RS 232
DOS20P	DOSP -2010	Platform Shaker with 20 mm linear motion
DOS20P	DSOP-2011	Platform Shaker with 20 mm linear motion along with RS 232
DOSUA-AB	Universal Attachment	Universal base with adjustable bars

Accessories

Part Name	Cat #	Details
Plate with slip resistance mat	DOSP - DA01	-
Funnel Attachment	DOSP - FA 01	-
Base Plate for Plastic clamp	DOSP -PPC	Base plate of 9 plastic clamps
Set of 9 x 250 ml Plastic clamp	DOSP -PS-250	-
Set of 9 x 500 ml Plastic clamp	DOSP -PS-500	-
Plastic clamp for 250 ml	DOSP -SPC-1	1 Plastic clamp - 250 ml
Plastic clamp for 500 ml	DOSP -BPC-1	1 Plastic clamp - 500 ml
Base plate for SS clamps 250 ml	DOSP -SSP-250	Base plate of SS clamp for 250 ml and smaller
Set of 9 x 250 ml SS clamp	DOSP -SS-250	-
Set of 13 x 150 ml SS clamp	DOSP -SS-150	-
Set of 13 x 100 ml SS clamp	DOSP -SS-150	-
Stainless steel clamp for 100 ml flask	DOSP -SC-100	Single clamp to hold 100 ml flask
Stainless steel clamp for 150 ml flask	DOSP -SC-150	Single clamp to hold 150 ml flask
Stainless steel clamp for 250 ml flask	DOSP -SC-250	Single clamp to hold 250 ml flask
Base plate for SS clamps 500 ml	DOSP -SSP-500	Base plate of SS clamp for 500 ml and higher
Set of 9 x 500 ml SS clamps	DOSP -SS-500	-
Set of 5 x 1000 ml SS clamps	DOSP -SS-1000	-
Stainless steel clamp for 500 ml flask	DOSP -SC-500	Single clamp to hold 500 ml flask
Stainless steel Clamp for 1000 ml flask	DOSP -SC-1000	Single clamp to hold 1000 ml flask

Contact our Sales Department on

+61 2 7233 3237

Sales@diaguru.com.au

www.diaguru.com.au