






Sharps container catalogue



	DESCRIPTION	CODE	CARTON QTY
	<p>FITTUBE™ Container Yellow 500ml Product summary: Compact high quality all-purpose injection moulded container Material: Polypropylene Nominal capacity: 450 ml (safe recommended usable space) Dimensions: 115mm high, 70mm x 85mm, 57mm lid Container shape: Rectangular Product standards: Australian Standards Accreditation - AS4031:1992 Manufacturing standards: ISO 9001:2015</p>	DS-201	Units per cartons: 133
	<p>FITTUBE™ Container Large - Yellow 750ml Product summary: Compact high quality all-purpose injection moulded Material: Polypropylene Nominal capacity: 675 ml (safe recommended usable space) Dimensions: 180mm high, 85mm x 70mm, 57mm lid Container shape: Rectangular Product standards: Australian Standards Accreditation - AS4031:1992 Manufacturing standards: ISO 9001:2015</p>	DS-202	Units per cartons: 84
	<p>FITTANK™ Container 1.4 Lt Snap Top with petal insert Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.2 litres (safe recommended usable space) Dimensions: 202.5mm (h), 97mm (l), and 91mm (w) Container shape: Rectangular Syringe Capacity: 40 x 1ml, 27 x 2ml, 14 x 5ml and 7 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-127	Units per cartons: 58
	<p>FITTANK™ Container 1.4 Lt Re-Sealable Top with petal insert and needle notch Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.2 litres (safe recommended usable space) Dimensions: 202.5mm (h), 97mm (l), and 91mm (w) Container shape: Rectangular Syringe Capacity: 40 x 1ml, 27 x 2ml, 14 x 5ml and 7 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-167	Units per cartons: 58
	<p>FITTANK™ Container 1.4 Lt Screw Top - no insert Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.2 litres (safe recommended usable space) Dimensions: 202.5mm (h), 97mm (l), and 91mm (w) Container shape: Rectangular Syringe Capacity: 40 x 1ml, 27 x 2ml, 14 x 5ml and 7 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-128	Units per cartons: 58

	DESCRIPTION	CODE	CARTON QTY
	<p>FITTANK™ Container 1.4 Lt Screw Top with petal insert and needle notch Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.2 litres (safe recommended usable space) Dimensions: 202.5mm (h), 97mm (l), and 91mm (w) Syringe Capacity: 40 x 1ml, 27 x 2ml, 14 x 5ml and 7 x 10ml Product standards: Australian Standards Accreditation - AS4031:1992 Manufacturing standards: ISO 9001:2015</p>	DS-131	Units per cartons: 58
	<p>FITTANK™ Container 1.8 Lt Snap Top with petal insert Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.5 litres (safe recommended usable space) Dimensions: 235.9mm (h), 97mm (l), and 91mm (w) Syringe Capacity: 55 x 1ml, 32 x 2ml, 16 x 5ml and 10 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-129	Units per cartons: 48
	<p>FITTANK™ Container 1.8 Lt Re-Sealable Top with petal insert and needle notch Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.5 litres (safe recommended usable space) Dimensions: 235.9mm (h), 97mm (l), and 91mm (w) Syringe Capacity: 55 x 1ml, 32 x 2ml, 16 x 5ml and 10 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-168	Units per cartons: 48
	<p>FITTANK™ Container 1.8 Lt Screw Top – no insert Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.5 litres (safe recommended usable space) Dimensions: 235.9mm (h), 97mm (l), and 91mm (w) Syringe Capacity: 55 x 1ml, 32 x 2ml, 16 x 5ml and 10 x 10ml Manufacturing standards: ISO 9001:2015</p>	DS-130	Units per cartons: 48
	<p>FITTANK™ Container 1.8 Lt Screw Top with petal insert and needle notch Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label Material: HD Polypropylene Nominal capacity: 1.5 litres (safe recommended usable space) Dimensions: 235.9mm (h), 97mm (l), and 91mm (w) Syringe Capacity: 55 x 1ml, 32 x 2ml, 16 x 5ml and 10 x 10ml Product standards: Australian Standards Accreditation - AS4031:1992 Manufacturing standards: ISO 9001:2015</p>	DS-132	Units per cartons: 48

	DESCRIPTION	CODE	CARTON QTY
	<p>FITTANK™ Container 3 Lt Resealable Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 123mm (h), 231mm (l), & 136mm (w)</p> <p>Syringe Capacity: Volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-163	Units per cartons: 20
	<p>FITTANK™ Container 3 Lt Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 123mm (h), 231mm (l), & 136mm (w)</p> <p>Syringe Capacity: Volume sharps disposal</p> <p>Product standards: Australian Standards Accreditation - AS4031:1992</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-188	Units per cartons: 20
	<p>FITTANK™ Container 5 Lt Resealable Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 275mm (h), 231mm (l), & 136mm (w)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-195	Units per cartons: 12
	<p>FITTANK™ Container 5 Lt Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 275mm (h), 231mm (l), & 136mm (w)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Product standards: Australian Standards Accreditation - AS4031:1992</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-196	Units per cartons: 12

	DESCRIPTION	CODE	CARTON QTY
	<p>FITTANK™ Container 8 Lt Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 231mm (l) x 136mm (w) x 305mm (h)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Product standards: Australian Standards Accreditation - AS4031:1992</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-197	Units per cartons: 9
	<p>FITTANK™ Container 8 Lt Re-Sealable Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 231mm (l) x 136mm (w) x 305mm (h)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-164	Units per cartons: 9
	<p>FITTANK™ Container 12 Lt Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 270mm (h) x 195mm (w)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Product standards: Australian Standards Accreditation - AS4031:1992</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-176	Units per cartons: 7
	<p>FITTANK™ Container 19 Lt Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 330mm (h) x 195mm (w)</p> <p>Syringe Capacity: Large volume sharps disposal</p> <p>Product standards: Australian Standards Accreditation - AS4031:1992</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-180	Units per cartons: 5

	DESCRIPTION	CODE	CARTON QTY
	<p>Sharps EasyCollect 12 Litre Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992 Material: Injection Moulded Polyolefin Normal Capacity: 8 Litres (safe recommended usable space) Capacity: 12 Litres (maximum usable space) Dimensions: 365mm (h) x 179mm (w) Syringe Capacity: Large volume sharps disposal Product standards: Australian Standards Accreditation - AS4031:1992 Manufacturing standards: ISO 9001:2015</p>	DS-227CTN	Units per cartons: 5
	<p>Sharps EasyCollect 12 Litre NO TRAY Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992 Material: Injection Moulded Polyolefin Normal Capacity: 10 Litres (safe recommended usable space) Capacity: 12 Litres (maximum usable space) Dimensions: 365mm (h) x 179mm (w) Syringe Capacity: Large volume sharps disposal Manufacturing standards: ISO 9001:2015</p>	DS-227-NT	Units per cartons: 5
	<p>Sharps EasyCollect 18 Litre Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992 Material: Injection Moulded Polyolefin Normal Capacity: 14 Litres (safe recommended usable space) Capacity: 18 Litres (maximum usable space) Dimensions: 505mm (h) x 179mm (w) Syringe Capacity: Large volume sharps disposal Manufacturing standards: ISO 9001:2015</p>	DS-228CTN	Units per cartons: 5
	<p>Sharps EasyCollect 18 Litre NO TRAY Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992 Material: Injection Moulded Polyolefin Normal Capacity: 16 Litres (safe recommended usable space) Capacity: 18 Litres (maximum usable space) Dimensions: 505mm (h) x 179mm (w) Syringe Capacity: Large volume sharps disposal Manufacturing standards: ISO 9001:2015</p>	DS-228-NT	Units per cartons: 5

	DESCRIPTION	CODE	CARTON QTY
--	-------------	------	------------



Sharps EasyCollect 25 Litre

Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992

Material: Injection Moulded Polyolefin

Normal Capacity: 21 Litres (safe recommended usable space)

Capacity: 25 Litres (maximum usable space)

Dimensions: 690mm (h) x 179mm (w)

Syringe Capacity: Large volume sharps disposal

Manufacturing standards: ISO 9001:2015

DS-229CTN Units per cartons: 3



Sharps EasyCollect 25 Litre NO TRAY

Product summary: Compact high quality all-purpose injection moulded disposal container & lid with standard warning label. Impact, crush and needle penetration exceed AS4031:1992

Material: Injection Moulded Polyolefin

Normal Capacity: 23 Litres (safe recommended usable space)

Capacity: 25 Litres (maximum usable space)

Dimensions: 690mm (h) x 179mm (w)

Syringe Capacity: Large volume sharps disposal

Manufacturing standards: ISO 9001:2015

DS-229-NT Units per cartons: 3



SHARPS BIN DUO - 30 Litre

Product summary: Compact high quality injection moulded container, the lid features a screw top opening which can be locked in place. This container provides storage for a wide range of sharps items.

Material: Polyolefin

Normal Capacity: 30 Litres

Capacity to fill line: 22.5 Litres

Dimensions: 390mm (h) x 300mm (w) x 400mm (L)

Syringe Capacity: Large volume sharps disposal

Product standards: UN Certification

Manufacturing standards: ISO 9001:2015

DS-216 Units per cartons: 5



SHARPS BIN DUO - 60 Litre

Product summary: Compact high quality injection moulded container, the lid features a screw top opening which can be locked in place. This container provides storage for a wide range of sharps items.

Material: Polyolefin

Normal Capacity: 60 Litres

Capacity to fill line: 52.5 Litres







Dimensions: 650mm (h) x 300mm (w) x 400mm (L)

Syringe Capacity: Large volume sharps disposal

Product standards: UN Certification

Manufacturing standards: ISO 9001:2015

DS-220 Units per cartons: 5

DESCRIPTION	CODE	CARTON QTY
 <p>FITTUBE™ Container 1.4 Lt Cytotoxic Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 97mm (l) x 91mm (w) x 202.5mm (h)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-199	Units per cartons: 58
 <p>FITTUBE™ Container 1.8 Litre Cytotoxic Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 97mm (l) x 91mm (w) x 235.9mm (h)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-200	Units per cartons: 48
 <p>FITTUBE™ Container 5 Litre Cytotoxic - Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 231mm (l) x 136mm (w) x 275mm (h)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-205	Units per cartons: 12
 <p>FITTUBE™ Container 8 Litre Cytotoxic - Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 231mm (l) x 136mm (w) x 305mm (h)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-203	Units per cartons: 9
 <p>FITTANK™ Container 12 Litre Cytotoxic Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 270mm (h) x 195mm (w)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-186	Units per cartons: 7
 <p>FITTANK™ Container 19 Litre Cytotoxic Screw Top with petal insert and needle notch</p> <p>Product summary: Compact high quality all-purpose blow moulded disposal container & injection moulded lid with standard warning label</p> <p>Material: HD Polypropylene</p> <p>Dimensions: 330mm (h) x 195mm (w)</p> <p>Syringe capacity: Large volume sharps disposal</p> <p>Manufacturing standards: ISO 9001:2015</p>	DS-185	Units per cartons: 5

ITEM CODE & DESCRIPTION	CARTON QTY
 <p>DS-B003: 500/750ml Bracket Suitable for: DS-201: Fitsafe Container 500ml DS-202: Fitsafe Container 750ml</p>	Units per carton:10
 <p>DS-B001: FITTANK™ Bracket – 1.4/1.8 Litre Suitable for: DS-167: FITTANK™ Container 1.4 Litre DS-131: FITTANK™ Container 1.4 Litre DS-168: FITTANK™ Container 1.8 Litre DS-132: FITTANK™ Container 1.8 Litre</p>	Units per carton:10
 <p>DS-B002A: Quick Release Bracket Suitable for: DS-163: FITTANK™ Container 3 Litre DS-188: FITTANK™ Container 3 Litre DS-195: FITTANK™ Container 5 Litre DS-196: FITTANK™ Container 5 Litre DS-164: FITTANK™ Container 8 Litre DS-197: FITTANK™ Container 8 Litre DS-176: FITTANK™ Container 12 Litre DS-180: FITTANK™ Container 19 Litre</p>	Units per carton:10
 <p>DS-B006: Hanger Bracket – 3 Litre Suitable for: DS-163: FITTANK™ Container 3 Litre DS-188: FITTANK™ Container 3 Litre</p>	Units per carton: 5
 <p>DS-B007: Hanger Bracket – 5 & 8 Litre Suitable for: DS-195: FITTANK™ Container 5 Litre DS-196: FITTANK™ Container 5 Litre DS-164: FITTANK™ Container 8 Litre DS-197: FITTANK™ Container 8 Litre</p>	Units per carton: 5
 <p>DS-B005: Bracket Basket – 30 Litre Suitable for: DS-216: Sharps Bin 30 Litre DUO</p>	Units per carton: 5
 <p>DS-B008: Bracket Easycollect Suitable for: DS-227: EasyCollect 12 Litre DS-227-NT: EasyCollect 12 Litre No Tray DS-228: EasyCollect 18 Litre DS-228-NT: EasyCollect 18 Litre No Tray DS-229: EasyCollect 25 Litre DS-229-NT: EasyCollect 25 Litre No Tray</p>	Units per carton: 10
<p>DS-B009: Trolley to Suit 30/50/60 Litre Suitable for: DS-216: Sharps Bin 30 Litre DUO DS-172: Sharps Bin 50 Litre DUO DS-220: Sharps Bin 60 Litre DUO</p>	Units per carton: 1
<p>DS-B010: Trolley EasyCollect Suitable for: DS-227: Sharp EasyCollect 12 DS-228: Sharp EasyCollect 18 DS-229: Sharp EasyCollect 25</p>	Units per carton: 1



WISHMED PTY LTD

14 Butterfield Street, Blacktown,
NSW 2148, AUSTRALIA

Telephone: 02 8678 6499

Email: sales@wishmed.com.au

www.wishmed.com.au